

# SENSITIVE SITE EXPLOITATION KIT

SPEC SHEET: TACTICAL ELECTRONICS

EOD TOOL  
KITS



## DESCRIPTION

The Sensitive Site Exploitation Kit is a small, lightweight but rugged system built with the SOF mission in mind. The kit provides the tools necessary to conduct numerous Sensitive Site Exploitation operations in an austere environment. Developed by operators with decades of SOF experience, this sampling system is all inclusive and organized into a small lightweight package. Two mini pelican cases are provided for you to customize your load out for each collection and also protect collected SSE material and evidence. A replenishment kit is also included, allowing you to easily restock your sampling kit.

The Sensitive Site Exploitation Kit provides the ability to exploit and safely sample:

- › DNA
- › Homemade Explosives (HME)
- › Chemical and Biological Warfare Agents (solid and liquid)
- › Drugs (various forms)
- › Safely remove and transport valuable electronics

## SPECIFICATIONS

<b>EXTERNAL BOX DIMENSIONS:</b>	19.0in x 14.25in x 5.75in
<b>REPLENISHMENT KIT DIMENSIONS:</b>	12.25in x 9.0in x 4.75in
<b>PELICAN CASES DIMENSIONS:</b>	9.88in x 5.59in x 2.63in

## FEATURES

- › Provides the ability to safely sample liquids, solids, and DNA for bulk or trace evidence in an austere environment
- › Modular, small, and lightweight kit designed to fit in any bag or case
- › Small outer case doubles as a hardened and airtight container to protect evidence, sample(s) or other SSE items like electronics
- › Small outer case provides a convenient and safe method of delivering hazardous samples to a laboratory, includes labeling system to provide key information to the lab technician
- › Sensitive Site Exploitation Kit provides a way to conduct light decontamination of samples or other materials, wanted in a biological or chemical environment
- › Sensitive Site Exploitation Kit provides items needed for a hasty decontamination station, such as a boot wash, drop cloths, and decon sprayer
- › Specialized sampling containers utilize varying barrier types for collecting liquids or solids that are reactive to standard bottles or vials
- › Sampling containers are specifically made for holding liquids that are caustic or acidic
- › Anti-static bags provided for handling materials that are sensitive to static electricity
- › Cost effective collection kit with the ability to perform multiple forensic collections or SSE operations
- › Curated by operators with decades of combined US/UK SOF experience



**CONTACT :**

TACTICAL ELECTRONICS  
P: 866 / 541 / 7996  
P: 918 / 249 / 8326  
F: 918 / 249 / 8328

**WEBSITE :** TACTICALELECTRONICS.COM

**DOCUMENT :**

OFFICIAL SPECIFICATIONS  
SHEET \*PROPRIETARY\*

**COPYRIGHT :**

©2018 TACTICAL ELECTRONICS  
AND MILITARY SUPPLY, LLC.

**PAGE NUMBER :**

PRODUCT SPECIFICATIONS  
AND INFORMATION ARE  
SUBJECT TO CHANGE  
WITHOUT NOTICE.

1

# SENSITIVE SITE EXPLOITATION KIT

SPEC SHEET: TACTICAL ELECTRONICS

EOD TOOL  
KITS



## EXTERNAL BOX

- › Outer Cardboard Box
- › Inner Cardboard Box (Replenishment Kit)
- › Decon spray bottle for bleach
- › (2) Foldable buckets, 10L
- › (4) Glass pipettes, 10 mL
- › 2 Pages of SSE Kit labels
- › Pelican Case #1
- › Pelican Case #2

## REPLENISHMENT KIT

- › (1) Jar of bleach tablets
- › (1) Roll of Gorilla tape
- › (8) Hype wipes
- › (8) Anti-static Ziploc bags, small
- › (8) Anti-static Ziploc bags, large
- › (8) Sets of Nitrile gloves, XL
- › (8) Pipettes, 7mL, Teflon, Karter Scientific
- › (2) Pipette pumps, Karter Scientific 10mL
- › (4) Syringes, 10cc
- › (8) IV extension tubes, luer lock, 21 inches
- › (4) Drop cloths, 4mil (3x2ft)
- › (1) Ph paper
- › (4) Centrifuge tubes, 50mL
- › (8) Wooden Q-tips with cap, 6in
- › (4) Blunt tipped needles, 22 gauge
- › (4) Glass 5mL bottles
- › (4) 30mL HDPE jars
- › (2) Luer locks
- › (8) Paper evidence bags
- › Parafilm® (8ft)

## PELICAN CASE #1

- › (1) Pelican™ 1060 case
- › (2) Hype wipes
- › (2) Anti-static Ziploc bags, small
- › (2) Anti-static Ziploc bags, large
- › (2) Sets of Nitrile gloves, XL
- › (1) Mini Sharpie pen
- › (2) Pipettes, 7mL, Teflon, Karter Scientific
- › (1) Syringe, 10cc
- › (2) IV extension tubes, luer lock, 21 inches
- › (2) Drop cloths, 4mil (3'x2')
- › (1) Telescopic alligator clip
- › (1) Set of Medical shears
- › (2) Blunt tipped needle 22 gauge
- › (2) 5mL bottle
- › (2) 30mL HDPE jar
- › (1) Lid Organizer
- › (1) Luer Locks
- › (2) Paper evidence bags
- › (2) Wooden Q-tips with cap, 6in
- › Parafilm® (1ft)



### CONTACT :

TACTICAL ELECTRONICS  
P: 866 / 541 / 7996  
P: 918 / 249 / 8326  
F: 918 / 249 / 8328

WEBSITE : TACTICALELECTRONICS.COM

### DOCUMENT :

OFFICIAL SPECIFICATIONS  
SHEET \*PROPRIETARY\*

### COPYRIGHT :

©2018 TACTICAL ELECTRONICS  
AND MILITARY SUPPLY, LLC.

### PAGE NUMBER :

PRODUCT SPECIFICATIONS  
AND INFORMATION ARE  
SUBJECT TO CHANGE  
WITHOUT NOTICE.

2

# SENSITIVE SITE EXPLOITATION KIT

SPEC SHEET: TACTICAL ELECTRONICS

EOD TOOL  
KITS



## PELICAN CASE #2

- › (1) Pelican™ 1060 case
- › (2) Hype wipes
- › (2) Anti-static Ziploc bags, small
- › (2) Anti-static Ziploc bags, large
- › (2) Sets of Nitrile gloves, XL
- › (1) Mini Sharpie pen
- › (2) Pipettes, 7mL, Teflon, Karter Scientific
- › (1) Syringe, 10cc
- › (2) IV extension tubes, luer lock, 21 inches
- › (2) Drop cloths, 4mil (3'x2')
- › (1) Telescopic alligator clip
- › (1) Set of Medical shears
- › (2) Blunt tipped needle 22 gauge
- › (2) 5mL bottle
- › (2) 30mL HDPE jar
- › (1) Lid Organizer
- › (1) Luer Locks
- › (2) Paper evidence bags
- › (2) Wooden Q-tips with cap, 6in
- › Parafilm® (1ft)



**CONTACT :**

TACTICAL ELECTRONICS  
P: 866 / 541 / 7996  
P: 918 / 249 / 8326  
F: 918 / 249 / 8328

**WEBSITE :** TACTICALELECTRONICS.COM

**DOCUMENT :**

OFFICIAL SPECIFICATIONS  
SHEET \*PROPRIETARY\*

**COPYRIGHT :**

©2018 TACTICAL ELECTRONICS  
AND MILITARY SUPPLY, LLC.

**PAGE NUMBER :**

PRODUCT SPECIFICATIONS  
AND INFORMATION ARE  
SUBJECT TO CHANGE  
WITHOUT NOTICE.

3